



diagnostics	87 - 100
blood sampling devices	88 - 93
urine testing	93 - 96
drug and allergy testing	97 - 99
pregnancy tests	97
nerve testing	100
urine and specimen containers	100

diagnostics



Cobas h232



Tests for diagnosis heart failure, DVT and myocardial infarction annotations.

- Rapid, on-the-spot decision support for treatment, referral, or discharge of cardiovascular patients
- Results available in 8 – 12 minutes
- Portable device which is suitable for use in the GP's surgery, one stop clinics or community hospitals
- No maintenance
- Easy to clean
- On board QC
- Results comparable to Roche laboratory method
- Quality assurance via patient identification and QC operator lock-out



NB. Delivery of BLD350 not included in price

BLD350	Cobas h232	£1361.72
BLD351	Cardiac D-Dimer tests, pack of 10	£82.37
BLD352	Cardiac D-Dimer control	£47.76
BLD353	Cardiac proBNP tests, pack of 10	£347.78
BLD354	Cardiac proBNP control	£86.94
BLD355	Cardiac T Quant tests, pack of 10	£140.05
BLD356	Cardiac Troponin-T control	£36.23
BLD357	Cardiac pipettes 150 UI, pack of 20	£16.05

Clinispin Horizon benchtop centrifuge



- Very simple to use
 - Quiet operation
 - Single speed 3,800 RPM/ 2,000 x g
 - Preset timer, variable 0-30 minute
 - Dimensions (hwxwd): 220 x 280 x 320mm
 - Supplied with 6-place swing-out rotor, 6 x 10ml adaptors & 6 x 5ml adaptors
 - Weight: 8.4kg
- NB. Delivery not included in price**

BLD521	£566.52
--------	---------

Safe-T-Pro Plus blood sampling device



- Single use, self destruct safety lancet
- Ideal for all multi-patient testing sites
- Choice of 3 depth settings
- Minimise the risk of needle stick injury
- Prevents cross-contamination between patients
- Quick and easy to use
- Box of 200

BLD124	Safe-T-Pro blood sampling device	£28.31
--------	----------------------------------	--------

Unistik 3, extra depth



- Suitable for tough skin and larger blood samples required for home cholesterol testing
- Single use safety lancet
- With Comfort Zone technology
- Designed to make blood sampling less painful
- Gauge: 21g (0.81mm)
- Depth: 2mm
- Pack of 100



DRE170

£9.60

Unistik 3, normal depth



- Unistik 3 normal lancets 23g/1.8mm
- For standard tests and normal skin types
- Safe and comfortable sampling for both yourself and others
- It is ideal if you are travelling, at work or at school or college
- Single use capillary blood sampling device, with the needlepoint hidden before use
- Unistik 3 makes use of Comfort Zone Technology, which masks the discomfort of lancing by targeting the sensation and perception of pain
- To avoid accidental injuries and the risk of cross infection, the needle retracts after use
- Gauge: 23g
- Depth: 1.8mm
- Medium flow
- Pack of 100



BLD6002

£9.60

Clearview D-dimer tests

The Clearview Simplify D-dimer test is a rapid test for the qualitative detection of D-dimer in human whole blood (either capillary or venous) and plasma. The test is an aid in the assessment and evaluation of patients with suspected disseminated intravascular coagulation (DIC), deep vein thrombosis (DVT) and pulmonary embolism (PE).

When used in conjunction with a validated pre-test probability score and diagnostic algorithm such as Wells, the Clearview Simplify D-dimer test can be used to safely rule-out VTE in low or moderate risk patients. It is estimated that only 20 – 30% of patients suspected actually have DVT.

Features

- Fast decisions for effective patient management
- Easily performed test without the need for instrumentation
- No specialised training required
- High negative prediction value (NPV) for the safe elimination of venous thromboembolism (VTE)



BLD360

£179.00

CoaguChek XS and CoaguChek XS Plus

CoaguChek XS

State of the art technology with the accuracy, safety and convenience essential for effective patient self testing.

- INR testing
- **Very high levels of accuracy** (ISI of 1.0) and precision (CV <5%)
- Meter provides the latest technology, quick and easy INR testing
- Heparin insensitivity
- Supplied with rechargeable battery pack, base unit and initial choice of lancing systems
- Compact, hand held design and battery operation
- Top and side dosing options with strip
- Individual therapeutic range
- Large graphical display and easy to follow operation
- 100 result memory

CoaguChek XS Plus

High specification meter for use in Primary care.

- Bringing anti coagulation services into GP Practices
- Full graphical touch screen for simple, fast handling
- Audit facility with patient and user identification
- Accurate results within one minute
- Download capabilities and integrated data management function
- Mains and battery power supply for in-clinic or home visit use
- 1000 result memory
- Level 4 funding available through National Enhanced Services
- Potential QOF Points available for XS Plus



A complete range for patients and practices



BLD460	XS Plus kit		£887.13
BLD465	XS patient meter	1	£473.48
BLD466	Test strips for XS Plus and XS, pack of 24	2	£77.51
BLD461	Test strips for XS Plus and XS, pack of 48		£151.54
BLD468	Control solution for XS Plus		£17.25
BLD124	Safe-T-Pro blood sampling device		£25.88

After purchasing your product please contact Roche Diagnostics to receive information about training and technical support services

Accutrend Plus cholesterol, glucose, triglycerides, and lactate testing system

- Simple to use and cost efficient
- Rapid results
- Accredited testing
- In addition to cholesterol testing, there are three other test parameters that can be run on the Accutrend Plus
- Cholesterol – quick, cost effective screening for cardiovascular risk
- Glucose – risk assessment for diabetes
- Triglycerides – high blood triglycerides are associated with an increased risk of developing coronary heart disease
- Lactate – used in sports therapy and in protocols such as sepsis and diabetic ketoacidosis
- Please note all items in KIT (BLD171) have a shelf life of 6 months



BLD171	Glucose, cholesterol, triglycerides & lactate kit	1	£276.99
BLD170	GC-cholesterol and glucose testing system	2	£212.81
BLD172	Triglycerides strips, pack of 25		£56.01
BLD173	Lactate strips, pack of 25	3	£43.44
BLD112	Cholesterol strips, pack of 25		£56.27
BLD056	Glucose strips, pack of 25		£8.91
BLD113	Cholesterol control solution	4	£11.03
BLD114	Glucose control solution		£10.48
BLD174	Triglycerides control solution		£18.98
BLD175	Lactate control solution		£18.98
BLD124	Safe-T-Pro blood sampling device		£25.88

CardioChek PA cholesterol and glucose profiler

The CardioChek PA is simple and convenient to use. It's the fastest, most cost-efficient and user-friendly method to accurately determine lipid values at point-of-care. In just 90 seconds, the CardioChek PA test system puts the specific patient data you need in hand. These attributes establish the CardioChek PA device as the standard in Point of Care cholesterol testing.



- The UK's number one choice for point-of-care lipid analysis
- Evaluates total cholesterol, HDL, LDL, triglycerides, glucose & ketones in various combinations
- Fully automated, one sample provides results for up to 4 different parameters
- The world's first portable direct LDL cholesterol test
- Results in 90 seconds
- Ideal for pharmacy, surgery, occupational health and high volume health screenings
- Accurate – CE marked/FDA cleared/CLIA waived/CRMLN certified
- No need to refrigerate strips
- Stores up to 200 results – date and time stamped

Single tests require 15µl pipette, double tests require a 30µl pipette, and triple/quad tests require a 40µl pipette

BLD580	CardioChek PA starter kit, includes 15 cholesterol/HDL/glucose	£799.00
BLD560	CardioChek PA monitor	£510.00
BLD562	15µl pipette, 15/pk	£7.91
BLD582	30µl pipette, 25/pk	£6.20
BLD562	40µl pipette, 16/pk	£7.91
BLD565	Glucose panels, 25/pk	£12.95
BLD566	Cholesterol panels, 25/pk	£69.95
BLD567	HDL panels, 25/pk	£72.36
BLD574	Direct LDL panels, 25/pk	£110.00
BLD568	Triglyceride panels, 25/pk	£75.00
BLD569	Ketone panels, 25/pk	£75.00
BLD573	Cholesterol/glucose combined panels, 25/pk	£66.95
BLD575	Cholesterol/HDL combined panels, 25/pk	£135.00
BLD576	Cholesterol/HDL/glucose combined panels, 15/pk	£98.40
BLD577	HDL/triglycerides/glucose combined panels, 15/pk	£82.00
BLD561	Cholesterol/HDL/LDL/triglyceride panels, 15/pk	£159.00
BLD570	Multichemistry control solution, 2 vials (high/low)	£49.95
BLD571	HDL cholesterol control solution, 2 vials (high/low)	£37.95

DCA Vantage analyser

The DCA Vantage analyser is the diabetes system that delivers the most proven, clinically relevant HbA1c and A:C ratio. It's simple to use and yet powerful enough to deliver lab-quality sensitivity and specificity.

- Monitor glucose levels with unparalleled precision, accuracy, and ease
- Detect even minor changes in HbA1c levels
- Microalbumin, creatinine, and A:C ratio results reported with one urine sample
- Automatic quantitative results; HbA1c results in 6 minutes; Microalbumin, creatinine and A:C ratio available in 7 minutes
- No sample or reagent preparation required
- Track patients more comprehensively and recommend the most effective treatment plan with new features like flexible patient result recording and patient trending graphs
- Store and organise up to 4000 results
- Onboard GFR calculator and other new features help to monitor and detect kidney disease earlier
- Data can be transferred easily to a PC using a USB flash drive
- Results can be reported onscreen, printed via integrated printer, external printer or uploaded to LIS/HIS
- Remote access for supervisors or POC coordinators
- Proven system with a 12-year track record
- Please contact our sales team for enquires about ordering reagent and control kits



SIEMENS

BLD180	DCA Vantage instrument	£3695.00
BLD181	DCA barcode reader	£260.96

STAT-Site M haemoglobin meter



- Test with children, infants and adults for conditions of anaemia, fluid retention, dehydration or recent haemorrhage
- Measures in mmol/l, µg/ml or mg/dl
- Self prompting graphics guides user through the simplest test procedure
- Results in seconds
- Recall last result
- Uses whole blood sample, finger prick
- No user calibration required
- Battery operated, batteries supplied
- **Test cards supplied in pack of 100 – not included with BLD550**


BLD550	Haemoglobin meter	£240.00
BLD551	Haemoglobin test cards, pack of 100	£125.00
BLD124	Safe-T-Pro blood sampling device	£28.31

Medi-Test urine analysis strips

Exceptionally high quality urine analysis strips for surgery or hospital use. Medi-Test has a higher sensitivity than many other brands and is ergonomically designed to make them particularly user friendly



- Increased sensitivity over leading brands, glucose sensitive to 1.1mmol/L
- Advanced colour stability provides reliable results for longer
- All results develop after 30 seconds for your convenience, no need to stagger readings
- User friendly, easy to read comparison chart makes anomalies stand out immediately
- Manufactured to exacting standards in Germany, a strip you can have confidence in



		Glucose	Protein 2	pH	Ketones	Blood	Nitrite	Bilirubin	Urobilinogen	Leucocytes	Specific gravity	Ascorbic acid	
DIS142	Glucose, pack of 50	•											£2.95
DIS143	Protein-2, pack of 50		•	•									£3.95
DIS144	Keton, pack of 50				•								£3.75
DIS14501	Combi-2, pack of 50	•	•										£7.95
DIS14502	Combi-2, pack of 100	•	•										£12.50
DIS146	Combi-5N, pack of 50	•	•	•		•	•					•	£9.95
DIS14701	Combi-5S, pack of 50	•	•	•	•	•							£9.95
DIS14702	Combi-5S, pack of 100	•	•	•	•	•							£17.95
DIS148	Combi 8, pack of 100	•	•	•	•	•	•			•	•		£18.95
DIS149	Combi 10SGL, pack of 100	•	•	•	•	•	•	•	•	•	•		£27.95

URYXXON Relax



Optional urine analysis strip reader for Medi-Test urine analysis strips.

- The URYXXON Relax makes urinalysis simpler and more reliable
- Fully customisable for SI, arbitrary (+++) or conventional units, patient and operator ID
- Provides accurate and dependable urine status results within 60 seconds
- Eliminating subjectivity of visual colour interpretation
- Built in printer for your records
- The comprehensive interface options and the optimized printouts minimize risks associated with manual transcriptions
- Reliable results can be obtained virtually immediately at the point of care
- Compatible with Medi-Test Combi-8 and Combi-10

NB. Delivery of DIS122, DIS124 or DIS117 not included in price

DIS122	Uryxxon Relax strip reader	£699.00
DIS124	Uryxxon paper rolls, pack of 5	£34.99
DIS117	Urinalysis control solution kit	£65.00
DIS132	Checkmate professional pregnancy tests, pack of 20	£29.00



Urisys 1100



Roche Diagnostics' Urisys 1100 analyser offers a user friendly, handy and time-saving way to perform fast, standardized urine test strip evaluations.

Results are printed at the integrated thermal printer, stored in the internal memory or easily transmitted to computer systems. The optional bar code reader allows fast and easy patient identification.

Suitable for rapid screening of diseases of the kidney and the urogenital tract, carbohydrate metabolism disorders (diabetes mellitus), liver diseases and haemolytic disorders.

- Compact and extremely easy to use
- Provides standardised results and virtually eliminates possible sources of visual evaluation errors
- Easy to operate – Incubation, evaluation and printout occur automatically, quickly and quietly
- The printer unit offers reprint functions
- Has a serial interface that allows easy connection to computer systems
- Other connectivity options include barcode/keyboard input of patient data
- Easy cleaning and changing of printer rolls

NB. Delivery of LAB150 not included in price

LAB150	Urisys 1100	£640.00
LAB153	Printer paper for Roche Urisys, pk 4 x 5 rolls	£19.00

Combur-Test urinalysis strips



Maximum accuracy, hygiene & value... minimum hassle & expense

- Resistant to vitamin C interference
- Hygienic to handle and easy to use
- All test zones read in 60 seconds
- Accurate, reliable results every time
- Excellent value for money



		Glucose	Protein	pH	Ketones	Blood	Nitrite	Bilirubin	Urobilinogen	Leucocytes	Specific gravity	Micro-albuminuria
DIS100	Combur-3 test pH, pack of 50 strips	•	•	•								£7.10
DIS101	Combur-4 test N, pack of 50 strips	•	•	•			•					£9.49
DIS103	Combur-7 test N, pack of 100 strips	•	•	•	•	•	•			•		£19.49
DIS118	Combur-9 test, pack of 50 strips	•	•	•	•	•	•	•	•	•		£20.50
DIS099	Combur-10 test, pack of 100 strips	•	•	•	•	•	•	•	•	•	•	£29.99

Clinitek Status+



- Printed record to help avoid transcription errors
- Ready for connectivity when you are – via a simple upgrade kit
- Comprehensive urine test menu including hCG pregnancy testing
- Detects common sample interferences
- Checks strip integrity for exposure to humidity and prevents reporting if exposure is detected
- Identifies type of strip automatically
- Provides simple, accurate results for higher standards in urinalysis
- Helps overcome possible human error inaccuracies in visually read results
- New level of precision and reliability
- Unprecedented convenience
- Reassuring proven performance
- Minimal training required and audit compliant
- Clinitek Microalbumin allows detection of kidney disease before it becomes irreversible
- On the spot semi-quantitative results in 1 minute
- Clinitek can be used with a wider range of multistix testing including: Multistix 10SG, Multistix 8SG, Multistix GP, Clinitek Microalbumin and Clinitest hCG, see Siemens diagnostic reagent strips below

SIEMENS

DIS108	Clinitek Status+	£749.00
DIS109	Clinitek printer paper, pack of 5 strips	£31.10
DIS110	Clinitek hCG, for Clinitek Status, pack of 25	£82.00
DIS111	Clinitek microalbumin, pack of 25 strips	£59.47

Diagnostic reagent strips

SIEMENS

- High quality reagent strips for urinalysis
- Can be used with the Siemens Clinitek Status+ urine analyser, shown above



	Glucose	Protein	pH	Ketones	Blood	Nitrate	Bilirubin	Urobilinogen	Leucocytes	Specific gravity
--	---------	---------	----	---------	-------	---------	-----------	--------------	------------	------------------

DIS112	Albustix, pack of 50 strips		•								£6.09
DIS113	Multistix GP, pack of 25 strips	•	•	•	•	•			•	•	£10.45
DIS114	Multistix 8SG, pack of 100 strips	•	•	•	•	•			•	•	£37.75
DIS107	Multistix 10SG, pack of 100 strips	•	•	•	•	•	•	•	•	•	£42.00
DIS139	Uristix, pack of 50 strips	•	•								£10.25
DIS116	Microalbustix, pack of 25 strips	•	•	•	•		•	•		•	£54.00
LAB227	Hema combistix euro, pack of 50 strips	•	•	•	•						£20.00
LAB229	Hemastix, pack of 50 strips						•				£18.54
LAB230	Labstix, pack of 100 strips	•	•	•	•					•	£33.00
LAB231	Labstix SG, pack of 100 strips	•	•	•	•					•	£39.00
LAB232	N-Labstix, pack of 100 strips	•	•	•	•	•					£38.00

One Step Plus urinalysis strip reading device

The HS OneStepPlus is a state-of-the-art urine test strip reader used to scan the results of the 'HS Urispec Plus 10SL' (DIS134). The measurements are performed under standardised conditions, and is designed for in-vitro diagnostic use.



- Quantitative results in 90 seconds which can be viewed on the touch screen, printed or transferred to a computer
- Stores 200 measurements which are automatically saved
- Patient identification can be entered via keyboard, barcode or via touch screen
- Password option for enhanced security
- Acoustic confirmation/warning of positive results
- Lot control feature to reduce risk of using expired test strips
- Adjustable level of sensitivity
- Auto start option
- Standby after 3 minutes of inactivity to save power
- Multilingual – choose from 9 languages
- Power options: Battery or mains
- In compliance with directives 89/336/EEC and 73/23/EEC



Plus 10SL

- Blood
- Protein
- Ketones
- pH
- Urobilinogen
- Nitrate
- Glucose
- Density
- Bilirubin
- Leucocytes

LAB200	One Step Plus urinalysis strip reading device	£595.00
LAB210	Barcode scanner	£199.00
LAB211	Printer paper, pack of 5	£20.55

DIS134	HS Urispec Plus 10SL, pack of 100 strips	£21.10
--------	--	--------

Urispec diagnostic test strips



Glucose	Protein	pH	Ketones	Blood	Nitrate	Bilirubin	Urobilinogen	Leucocytes	Specific gravity	Density
---------	---------	----	---------	-------	---------	-----------	--------------	------------	------------------	---------

DIS136	Plus glucose, pack of 100 strips	•									£7.70
DIS133	Plus 9+Leuco, pack of 100 strips	•	•	•	•	•	•	•	•	•	£18.99
DIS134	Plus 10SL, pack of 100 strips	•	•	•	•	•	•	•	•	•	£19.00
DIS135	Plus 11 SYS, pack of 100 strips	•	•	•	•	•	•	•	•	•	£22.15
DIS137	Plus 3, pack of 100 strips	•	•							•	£9.25
DIS138	Plus 5+Leuco, pack of 100 strips	•	•		•	•	•	•			£11.85





HS One Step pregnancy tests

One Step HCG urine test is a rapid pregnancy test, which you can easily carry out yourself. It detects the presence of human chorionic gonadotropin (hCG), which appears in urine very early during pregnancy.

- For external use only
- For professional use only
- Available as:
strip tests (pack of 50),
cassette tests (pack of 25) and
midstream tests (pack of 25)



LAB212	Urine test strip, pack of 50	1	£18.00
LAB213	Urine test cassette, pack of 25	2	£12.00
LAB214	Midstream urine test, pack of 25	3	£40.00

Biocard coeliac test



- For the diagnosis of coeliac disease
- Requires a small finger prick sample
- Results in 5 minutes
- Detects high levels of IgA anti-transglutaminase antibodies
- Single use test kit
- Accurate and reliable

BLD586	£19.99
--------	--------

GlutenTox professional gluten detection kit

Designed specifically for sterile environments and commercial kitchens, GlutenTox pro is a user-friendly gluten detection kit for all surfaces as well as foods and drinks. Comes in a practical and simple test strip format.

- Quick, user-friendly test to detect gluten in foods, drinks, cosmetics etc.
- Empowers and informs coeliac sufferers about the foods they consume
- The only test that specifically recognizes the toxic 33-mer fraction of gluten, giving less false positives
- Sensitive enough to distinguish between different strains of the same grain and different brands of the same foodstuff
- Adjustable sensitivity: test for 20ppm or 5ppm
- No special equipment required; everything contained in the kit
- Safe, non-toxic, and easy to use
- Quick response, results in less than 20 minutes



BLD587	£53.32
--------	--------

Fastcheck allergy test

Fastcheck is an excellent first step in diagnosing IgE mediated allergies with the minimum patient discomfort. This single use, fully disposable enzyme-linked immunoassay (ELISA) test enables allergy testing in almost any location.

- IgE antibody blood test for more than 24 different allergens, either 12 airborne allergens or 12 food allergens
- Covers 90% of the UK's most common allergens
- Safer and more comfortable than conventional allergy tests – patient has no contact with test allergens virtually eliminating risk of anaphylaxis
- Results visible to the naked eye in just 30 minutes
- Requires just 2 to 3 drops of capillary blood taken from fingertip
- Developed from established laboratory technology resulting in excellent performance
- A great first step in diagnosing and treating allergies
- Specific IgE blood testing is recommended by N.I.C.E. (GC116)
- Not suitable for children under 5 years old



Fastcheck inhaled

Checks for:

- House dust mite
- Cat epithelia
- Sting nettle pollen
- Hazel pollen
- Birch pollen
- Cladosporium herbarum
- Grass pollen mix
- Latex
- Aspergillus fumigatus
- Mugwort pollen
- Dog epithelia
- Alternaria alternata



BLD584	Fastcheck allergy test x 12	£29.99
--------	-----------------------------	--------

Cerascreen food intolerance test

Cerascreen is an excellent first step in diagnosing IgG food intolerances with the minimum patient discomfort. This single use, fully disposable enzyme-linked immunoassay (ELISA) test enables food intolerance testing in almost any location.

- IgG antibody blood test enables identification food intolerance after common symptoms, including bloating, nausea, diarrhoea and stomach pain
- Tests for the UK's 12 most common intolerances including gluten, milk and peanut
- More comfortable than conventional intolerance testing
- Results visible to the naked eye in just over 30 minutes
- Requires just 1 drop of capillary blood taken from fingertip
- Developed from established laboratory technology resulting in excellent performance
- A great first step in diagnosing and treating dietary intolerance
- Not suitable for children under 5 years old



Nutrition intolerances:

- Hen egg
- Potato
- Wheat flour
- Rice
- Peanut
- Fruits (mix)
- Meat (mix)
- Soy bean
- Gluten free grains
- Grains containing gluten
- Cow's milk & other milk



BLD585	Cerascreen food intolerance test	£45.00
--------	----------------------------------	--------



HS One Step faecal occult blood test

The One Step faecal occult blood (FOB) test strip is a simple direct binding immunoassay for the rapid, qualitative detection of human occult blood in faeces. The assay is designed to detect lower levels of colorectal bleeding than other biochemical methods.

NEW

- For in vitro diagnostic use only
- For professional use only
- Available as strip tests (pack of 50) and cassette tests (pack of 25)



LAB215	Test strip, pack of 50	1	£48.00
LAB216	Test cassette, pack of 25		£33.00

HS One Step multi-5 drugs of abuse panel

One Step multi-drug urine test panel offers any combination from of abuse tests for 5 different drugs: amphetamine (AMP), benzodiazepines (BZO), cocaine (COC), marijuana (THC), opiate (OPI 2000).

This is a rapid one step test for the qualitative detection of drug use and their principal metabolites in human urine at a specified cut off level.

- For professional use only
- For in vitro diagnostic use

NEW



LAB217	Pack of 20		£61.00
--------	------------	--	--------

HS One Step strep A cassette test

HS One Step strep A swab test is a rapid, one step lateral flow immunoassay for the qualitative detection of the Group A streptococcal antigen from throat swab specimens to aid in the early diagnosis of group A streptococcal infection.

- For in vitro use only

NEW



LAB218	Pack of 25		£36.00
--------	------------	--	--------

diagnostics

Neuropen

Neurotips™




Neuropen

- Test for peripheral neuropathy in diabetic patients
- Perfect dual purpose, pocket sized device
- Calibrated neurotips test assesses reduced sensation to sharpness/pain in small nerve fibres and touch/pressure sensation in large nerve fibres
- Each Neuropen comes complete with 1 x 10g monofilament, 1 Neurotip and user instructions

Neurotips

- Sterile, single-use neurological examination pins
- Calibrated force of 40g may be safely exerted onto the skin to identify patients with a decrease in sensation to sharpness/pain in small nerve fibres
- Neurotips are supplied in boxes of 100 containing 50 red and 50 white which provide a contrast when used for peripheral vision testing

PER021	Neuropen	£16.28
PER022	Calibrated monofilament, 10g x 5	£14.36
PER023	Neurotips, box of 100	£9.61



Bailey's retractable calibrated monofilament

The Bailey retractable monofilaments are diagnostic devices that consist of a monofilament mounted in a hand unit for protection when not in use.

The monofilaments are primarily used for testing the feet of diabetic patients for loss of protective sensation associated with peripheral neuropathy. If a patient cannot feel 10g, they should be considered at risk for developing pressure or neuropathic ulcers and should be referred to a specialist for further assessment.

- Supplied with full user instructions
- Easy to clean, ergonomic design for ease of use
- Convenient pocket clip



PER024	10g	£ 14.99
--------	-----	---------

Urine cups

- 125cc graduated unbreakable cup with frosted writing field

DIS615	Snap-on red cap, non-sterile, pack of 100	£17.50
DIS616	Screw-on green cap, non-sterile, pack of 100	£17.50

